




INSTITUTE OF CHEMICAL TECHNOLOGY

Department of Food Chemistry and Analysis

Prague, Czech Republic, jana.hajslova@vscht.cz

The Use of DART MS for Characterizing Food

Jana Hajslova



10th joint CSL / JIFSAN symposium on Food Safety and Nutrition, **Methods and Systems for Tracking, Tracing, and Verifying Foods**



May 13 – 15, 2009, Greenbelt Marriot Hotel, MD, USA

OVERVIEW

➡ **Setting the scene:**

- A lot of requirements for fast food analysis

➡ **Ambient mass spectrometry:**

- What is the principle and scope
- DART TOF MS - a challenging alternative

➡ **Examples of DART -TOF MS use:**

- Tracking / authentication
- Food safety control - traceability



A lot of analytical work on food quality / safety


- Proteins, lipids, saccharides
- Vitamins, minerals
- Flavours and odours
- Food additives, supplements
- Residues and contaminants
- Processing and packaging contaminants
- Mycotoxins, marine and plant toxins
- Genetically modified organisms (GMO's)
- Allergens
- Nanoparticles



- Compliance with legislation
- Label declarations
- Authenticity, traceability, fraud

Setting the scene:

**DART- TOF MS JOINS SPECTROSCOPIC
TRACKING TOOLS FAMILY**



Outline of the TRACE study



Trappist beer is a Belgium (very nice but rather expensive) speciality produced by monks

400 beer samples involved in study:



Rochefort
other Trappist
non-Trappist



► The aim of study:

DISTINGUISH ROCHEFORT 8° FROM OTHER TRAPPIST AND NON-TRAPPIST BEERS

CHEMOMETRIC ANALYSIS OF DART-TOF MS DATA


tools ► Linear Discriminant Analysis (**LDA**)
Artificial Neural Networks (**ANN**)

Trappist Rochefort 6°, 8°, 10° beers vs. rest of beers

LDA model

Recognition ability: **98%**

Prediction ability: **97%**




ANN model

Recognition ability: **100%**

Prediction ability: **98%**


m/z

DART TOF MS enables fast tracking of beer origin



No problem to distinguish these to Budweisers within few minutes!


Beer consumption (l)






AMBIENT MASS SPECTROMETRY

...opens the doors to many challenging applications in various areas including food / feed control



- **Reduced / minimal sample prep**
- **No chromatographic separation**



TRAC, Volume 27, Issue 4, April 2008, Pages 284-290

Recent review

Ambient desorption ionization mass spectrometry

Andre Venter, Marcela Nefliu, R. Graham Cooks

The ambient ionization methods retain the signature advantages of MS




➔ **speed**

➔ **chemical specificity**

➔ **low detection limits**


AMBIENT DESORPTION IONIZATION METHODS





a Techniques where **ESI** mechanisms are mainly responsible for ionization.

b Methods where **chemical ionization** is responsible for ionization
(photoionization - *PI*, ion evaporation - *IE* and electrical discharge)

DART ion source: → → → Direct Analysis in Real Time



Autosampler



sample acquisition

DART Ionization: ANALYTES

Negative ions

- Direct (M^-)
 - ionic compounds, some electrophiles
- Proton abstraction $[M-H]^-$
 - acidic compounds, nitroaromatics
- Adduct formation $[M+X]^-$
 - Unstable nitro compounds, some halocarbons

Positive ions

- Direct (M^+ , $M^{\cdot+}$)
 - ionic compounds, low-IP organics
- Proton transfer $[M+H]^+$
 - Bases, alkenes, small alcohols, ethers, ketones, aldehydes
 - H/D exchange
- Other adducts $[M+Z]^+$
 - Polar compounds, ethers, ketones, acids, peroxides

When interfaced to a high-resolution TOFMS → measurement of exact masses and accurate isotopic abundances → **identification of unknowns**

Case study # 1:

HONEY



- *What is the origin?*
- *Is it authentic?*



STRATEGIES IN ASSESSMENT OF HONEY QUALITY AND AUTHENTICITY

Target analysis of physico-chemical parameters / components / enzymes


Molecular biology, PCR (pollen)

NON-TARGET PROFILING

- ➔ Near infrared spectroscopy
- ➔ Fourier-transform mid-infrared spectroscopy
- ➔ Fourier-transform Raman spectroscopy
- ➔ Nuclear magnetic resonance spectroscopy
- ➔ Gas chromatography-mass spectrometry
- ➔ Direct Injection mass spectrometry

AMBIENT MASS SPECTROMETRY

Can direct be honey analysed directly?



NO!

High intensity of non-characteristic ions - mostly degradation products of sugars (carameliation)

DART: positive ion mode, DART temperature: 250°C

Honey + water (2 g + 2 mL)




PRE-CONCENTRATION USING SPME

➔ *isolation of volatile compounds in the presence of (abundant) non-target compounds (e.g. sugars)*

or

➔ *absorbtion also less volatile (more polar) fraction by immersion of a fibre into sample solution?*



HEADSPACE: honey + water (2 g + 2 mL)

DIRECT IMMERSION: honey + water (2 g + 6 mL)


DART TOF MS

Head-space SPME profile

Fiber:
DVB/CAR/PDMS 50/30 μm
(divinylbenzene/Carboxen/
polydimethylsiloxane)

DART: positive ion mode
DART temperature: 250°C

*lower MW compounds
extracted*






DART TOF MS direct immersion SPME profile

Fiber: DVB/CAR/PDMS 50/30 μm
(divinylbenzene/Carboxen/
polydimethylsiloxane)

DART: positive ion mode
DART temperature: 250°C


*Larger range of
compounds extracted*

Honey A




MS[1];1.045..1.106;-1.0*MS[1];0.759..0.779; /ESI+ / 20080328_SMPL1
 $\times 10^3$ Area (90742)


Honey B



in 10 min
AUTHENTIC HONEY
SAMPLES




Lime





Acacia



PHENOLICS




Rape



Sunflower

CLASSIFICATION USING PCA (20 marker masses selected)



Case study # 2:

OLIVE OIL



- *What is its origin?*
- *Is it authentic?*



Examined samples

various quality grade and botanical origin:

- Extra virgin olive oil (EVOO)
- Olive oil (OO)
- Olive pomace oil (OPO)
- Hazelnut oil (HO)



EXTRA VIRGIN OLIVE OIL

The oil obtained from the fruit of the olive tree solely by mechanical or other physical means under the conditions, particularly thermal conditions, that do not lead to alternations in the oil, and which has not undergone any treatment other than washing, decantation, centrifugation and filtration.

SAMPLE PREPARATION

- **TAGs analysis:** oil dilution with toluene (1:50, v/v)
- **Polar compounds analysis:** 2 min shaking of oil with MeOH–H₂O mixture (80:20, v/v)

DART–TOFMS method

- IONZATION MODE: positive
- ANALYSIS TIME: **1 min**
- GAS BEAM TEMPERATURE
 - (i) 350°C,
 - (ii) 220°C
- For TAGs analysis ammonia solution was employed as dopant




DART-TOFMS [+] mass spectra of **DILUTED OILS**

Extra virgin olive oil


FRAGMENT IONS

FRAGMENT IONS


[TAG+NH₄]⁺




Hazelnut oil







Olive pomace oil



Olive oil




DART-TOFMS [+] mass spectra of POLAR EXTRACTS



GROUPING ANALYSIS USING LINEAR DISCRIMINANT ANALYSIS (LDA)

■ Based on profiles of **TAGs**




- △ EVOO
- HO
- ◇ OPO
- OO
- ◆ MIX

- EVOO, HO, OO, OPO, MIX → *objects*
- TAGs: markers → 11 masses → *variables*

Prediction ability 100% for the EVOO/HO mixtures in the range **50:50 – 85:15 (v/v)**

GROUPING ANALYSIS USING LDA

Based on profiles of **POLAR COMPOUNDS**



△ EVOO
□ HO
◇ OPO
○ OO
◆ MIX

- EVOO, HO, OO, OPO, MIX → *objects*
- Polar compounds: markers → 12 masses → *variables*

Software statistiXL 1.8

Prediction ability 100% for the EVOO/HO mixtures in the range **50:50 – 94:6 (v/v)**

Note: Prediction ability was obtained on the basis of leave-one-out cross validation (LOOCV)

Accepted Manuscript

Title: Ambient mass spectrometry employing direct analysis in real time (DART) ion source for olive oil quality and authenticity assessment

Authors: Lukas Vaclavik, Tomas Cajka, Vojtech Hrbek, Jana Hajslova

PII: S0003-2670(09)00602-3
DOI: doi:10.1016/j.aca.2009.04.043
Reference: ACA 229908



Case study # 3:

ANIMAL FATS



- *Pork mixed with beef?*
- *Tallow added to lard?*



Classic approach: Fatty acids profile

Typical fatty acid composition of LARD and TALLOW

Transesterification

GC-FID



FATTY ACID	TALLOW (%)	LARD (%)
Lauric	1,0	< 1,0
Myristic	1,4-7,8	0,5-2,5
Palmitic	17,0-37,0	20,0-32,0
Palmitoleic	0,7-8,8	1,7-5,0
Stearic	6,0-40,0	5,0-24,0
Oleic	26,0-50,0	35,0-62,0
Linolic	0,5-5,0	3,0-16,0
Linolenic	< 2,5	< 1,5
Arachidonic	< 0,5	< 1.0
Eikosenoic	< 0,5	< 1,0


Lard and tallow differ
in TAGs composition ►

**Target markers for
DART**




TAG	TALLOW (%)	LARD (%)
LLO	0,0	2,6
PoOL	0,0	1,6
LLP	0,0	2,2
MOL	0,0	1,1
OOL	0,0	5,0
PoOO	2,1	2,0
POL	3,8	10,0
PoPO	5,9	2,9
MOP	5,1	0,0
PLP	0,0	2,6
OOO	3,4	5,7
POO	23,0	20,8
PLS	0,0	6,1
POP	10,8	8,0
MPS	3,7	0,0
SOO	11,4	4,6
POS	14,9	15,5
PPS	5,7	3,6
SOS	6,6	2,1
PSS	3,7	3,7

DART [+], m/z 200-925




LARD



TALLOW

CHEMOMETRIC ANALYSIS - LDA



- Lard, Tallow, MIX → *objects*
- TAGs: markers → 24 masses → *variables*


Software statistiXL 1.8

Model prediction ability
100% for the L / T mixtures
in the range **50:50 – 95:5 (v/v)**

Note: Prediction ability was obtained on the basis of leave-one-out cross validation (LOOCV)


DART [+], m/z 100-1100

0.297..0.382;-1.0*MS[1];0.114..0.253; /ESI+ / 20081030_maslo_OK




AUTHENTIC BUTTER

0.267..0.318;-1.0*MS[1];0.174..0.243; /ESI+ / 20081030_mleko_TUK



MILK FAT (control)


0.143..0.190;-1.0*MS[1];0.038..0.108; /ESI+ / 20081030_maslo_FAKE



COUNTERFEIT BUTTER


Case study # 3:

MELAMINE




- *How fast can be the baby food formula checked ?*
- *Can be the melamine content determined accurately?*





► Addition of melamine with high content of nitrogen in molecule increases apparent protein content

**Crystals with cyanuric acid
→ kidney damage, renal failure**




Courtesy of Dr. Dalen W. Agnew/Michigan State University



World Health Organization

TDI 0.5 mg/kg body weight



Crystals are visible in this histologic section of kidney from a cat with a history of eating some of the pet food that manufacturers recalled because of adulteration with melamine and co-contaminants. Toxicologic analysis of the kidney found cyanuric acid, melamine, and ammeline. Melamine and cyanuric acid can combine to form crystals.

Current analytical tools employed for determination of melamine in foods / feeds

LC-MS/MS

extraction
SPE clean-up
microfiltration/centrifugation
chromatogr. separation (HILIC)

LOQs: 5 – 250 ppb

GC-MS(/MS)

extraction
SPE clean-up
derivatization
chromatogr. separation

LOQs: 100 - 2000 ppb

ELISA (AgraQuant®,
Romer Labs)
extraction
incubation




LOQs: 0.1 – 250 ppm

MALDI-TOFMS
extraction
matrix preparation

SCREENING


**RAMAN
SPECTROSCOPY**
extraction
LODs: 0.05–0.1%

DART -TOF MS positive mass spectra: solvent standard and $^{13}\text{C}_3$ -melamine, 0.1 $\mu\text{g}/\text{ml}$



Can the sample prep be omitted?

**DART MS+ SPECTRUM:
MILK SPIKED WITH MELAMINE (5ppm)**






- ▶ extraction/clean-up needed
- ▶ instrument sensitivity should be improved


DART-TOFMS PERFORMANCE IMPROVEMENT by VAPUR™ API INTERFACE (IonSense, USA)

Principle: ions formed during DART ionization process are collected and transferred through the ceramic tube into the MS inlet



→ **effective ion transport** → **BETTER SENSITIVITY AND REPEATABILITY**



SIMPLIFIED SAMPLE PREPARATION PROCEDURE




DART TOF MS analysis of dried milk




Ion chromatogram (m/z 127.07): dried milk, spike 1.5 ppm




ass Chrom[1];127.01400..127.03200; / ESI+ / 20090124_SMPL1



INJECTION (1)




INTERFERENCE



QUANTIFICATION: isotope dilution technique

CALIBRATION PLOT



SOLVENT STANDARD

Matrix suppression

**MATRIX-MATCHED
STANDARD (dried milk)**

DART TOF MS method performance characteristics

PARAMETER	WITHOUT VAPUR	WITH VAPUR
LOD (S/N 3)	700 - 1000 ppm	100 – 150 ppm
LOQ (S/N 9)	1500 – 2800 ppm	250 – 280 ppm
REPEATABILITY (0.5 ppm)	10%	3 – 5%
RECOVERY (1.0 ppm)		98%

**Improved parameters by
VAPUR**



Real life samples: INTERLABORATORY COMPARISON


- Results of DART–TOFMS analyses were compared with those obtained by LC-MS/MS method (results provided by Eurofins)

SAMPLE	DART - TOFMS	RSD	LC-MS/MS
Condensed milk	4.04 ppm	3.1%	4.00 ppm
Dried milk (1)	2.33 ppm	3.5%	2.40 ppm
Dried milk (2)	0.51 ppm	4.2%	0.57 ppm




▶ CYANURIC ACID can be determined in the same sample

DART negative mass spectrum: spike 2 ppm




Case study # 4:

ACRYLAMIDE



- *Is it possible to determine acrylamide precursors?*
- *Can be predicted acrylamide formation extent in potato chips?*




MODEL SYSTEMS

- Solutions of glucose or fructose with asparagine in methanol:water (1:1, v/v)
- Quinine was added to each model mixture prior analysis (10µg/ml) as internal standard (I.S.)

Model	Asparagine (mg/ml)	Glucose (mg/ml)	Fructose (mg/ml)
M1	0.3	0.05	-
M2	0.3	0.15	-
M3	0.3	0.30	-
M4	0.3	-	0.05
M5	0.3	-	0.15
M6	0.3	-	0.30


MODEL M1 – M3, DART [+], 250°C

INCREASING GLUCOSE CONTENT




MODEL M4 – M6, DART [+], 250°C

INCREASING FRUCTOSE CONTENT



ACRYLAMIDE FORMATION IN MODEL SYSTEMS




→ Intensity of acrylamide was normalized to I.S. intensity (n=6)...

1 min. per sample !

Model	Asparagine (mg/ml)	Glucose (mg/ml)	Fructose (mg/ml)
M1	0.3	0.05	0.00
M2	0.3	0.15	0.00
M3	0.3	0.30	0.00
M4	0.3	0.00	0.05
M5	0.3	0.00	0.15
M6	0.3	0.00	0.30




**OK, it works in model , but
what about real-life
samples?**




POTENTIAL OF POTATOES TO YIELD ACRYLAMIDE

- ▶ Potatoes ($n=14$), various varieties
- ▶ Glucose, fructose, sucrose and asparagine were determined in methanol:water extracts using **HPLC–RID** and **HPLC–FLD**, respectively
- ▶ Formation of acrylamide during desorption/ionization process was monitored with **DART–TOFMS**





POTATO TUBER EXTRACT, DART [+], 250°C




METABOLOME !!!

POTATO EXTRACT, DART [+], 250°C




ACRYLAMIDE [M+H]⁺


INTERFERENCE




Normalized signal of acrylamide (m/z 72.045) versus glucose and fructose content in potatoes



**OK, but confirmation of
observed trends is needed
using established method**



Correlation of DART normalized acrylamide signal with acrylamide content in fried potato crisps determined by LC-MS/MS



SUMMARY

DART is primarily **small-molecules** analysis technique (with some exceptions) allowing rapid (real-time) **METABOLOMIC PROFILING (FINGERPRINTING)**

sample handling minimised / eliminated → workload decreased → analyses throughput increased

When DART coupled with HR TOF MS, identification of unknowns possible

**DART IS A REAL CHALLENGE IN
RAPID FOOD / FEED ANALYSIS**



Interested in recent food analysis innovations...?



4th International Symposium on

RECENT ADVANCES IN FOOD ANALYSIS

4–6 November 2009

Prague, Czech Republic

www.rafa2009.eu



- Residues and contaminants
- Authenticity, traceability, fraud
- Flavours and odours
- Processing and packaging contaminants
- Mycotoxins, marine and plant toxins
- Allergens
- Genetically modified organisms (GMO's)
- Nanoparticles
- Novel foods, nutritional supplements, organic food



SCIENTIFIC COMMITTEE (tentative list)

Prof. Jana Hajslova	<i>Institute of Chemical Technology, Prague, CZ (chair)</i>
Prof. Michel Nielen	<i>RIKILT-Institute of Food Safety, Wageningen, NL (co-chair)</i>
Prof. John Gilbert	<i>Central Science Laboratory, York, UK</i>
Dr. Samuel Godefroy	<i>Health Canada, Ottawa, Canada</i>
Prof. Hans-Gerd Janssen	<i>Unilever Research and Development, Vlaardingen, NL</i>
Prof. Rudolf Krska	<i>University for Agriculture and Applied Life Sciences, Tulln, A</i>
Dr. Steve Lehotay	<i>United States Department of Agriculture, Wyndmoor, USA</i>
Prof. Peter Schieberle	<i>Technical University of Munich, Garching, D</i>
Dr. Richard Stadler	<i>Nestlé Research Centre, Orbe, CH</i>
Prof. Franz Ulberth	<i>JRC, Institute for Reference Materials and Measurements, Geel, B</i>